

Development and validation of a point-of-care, urine assay to measure adherence to PrEP and ART

Giffin Daughtridge¹, Shane Hebel¹, Mark Fischl¹, Jumana Hashim¹, Elijah Kahn-Woods¹, Keith Kardos¹
1 – UrSure, Inc

Contact: giffin@ursureinc.com

Background

- PrEP and ART are highly effective at preventing new infections and suppressing HIV, but only when taken consistently¹
- Adherence to PrEP and ART is sub-optimal and current monitoring methods are inadequate to reach global scale²
- **UrSure developed a qualitative, visually-read, point of care urine test which can measure adherence to tenofovir (TFV)**, a metabolite of tenofovir disoproxil fumarate (TDF) and tenofovir alafenamide (TAF), prodrugs of most ART and all PrEP regimens
- The test uses lateral flow immunoassay, and has been validated with clinical samples against a LC-MS/MS machine

Methods

- We synthesized a novel derivative of TFV, conjugated that to BSA, and use that to develop and screen for an anti-TFV antibody
- Our anti-TFV monoclonal antibody yielded 100% sensitivity, and 97% specificity, with minimal cross-reactivity
- LFIA strip was optimized to the device’s cut off, based off internal analysis
- TFV concentrations above or below the cutoff were tested on the POC LFIA for 199 urine samples and validated against a LC-MS/MS test
- Storage and temperature stability tests were carried out to ensure integrity of the LFIA at +25% and -50% of the cutoff

Results

UrSure developed a POC LFIA prototype (Figure 1), which tests for TFV concentrations in urine at a set cutoff. **The prototype has 100% sensitivity and 100% specificity** against a urine LC-MS/MS TFV test. (Figure 2) We saw no changes in the performance of the LFIA with storage at room temperature, 45°C, and 55°C for up to 21 days. (Figure 3)

Figure 1. UrSure’s Rapid Tenofovir Test Prototype, ready for investigational use





Figure 2. Sensitivity, Specificity, of UrSure LFIA prototype

	LC-MS (+)	LC-MS (-)
LFIA (+)	37	0
LFIA (-)	0	162

Figure 3. Stability results in time and temperature study



Conclusions

- We have developed the first-ever urine TFV LFIA prototype with 100% sensitivity and specificity, available for investigational use.
- **UrSure’s POC TFV Test can inexpensively facilitate real-time monitoring of adherence to PrEP and TFV-based ART, helping providers to optimally allocate adherence resources**
- Due to the routine collection of urine in visits, the test is non-invasive, and seamlessly fits into clinical workflow.
- Enhanced adherence support can prevent seroconversions and improve outcomes for PrEP and ART patients

Next Steps

- Our prototype is currently going through further validation testing to ensure stability, precision, reproducibility and integrity in the presence of potentially interfering substances.
- UrSure’s Rapid Tenofovir Test will be used in over ten planned field settings for research studies around the world.

Acknowledgements

- Diagnostic Consulting Network for our LFIA development.
- NIH for SBIR Phase I (R43 MH114864) and II funding (R44 MH114864)

References

1. Kearney BP, Flaherty JF, Shah J. Tenofovir disoproxil fumarate: clinical pharmacology and pharmacokinetics. Clin Pharmacokinet. 2004;43(9):595–612
2. Iacob, S., et al. (2017). "Improving the Adherence to Antiretroviral Therapy, a Difficult but Essential Task for a Successful HIV Treatment—Clinical Points of View and Practical Considerations."