

Development of a Point of Care Urine Tenofovir Test – A PrEP and ART Adherence tool



Giffin Daughtridge, MD,MPA¹, Graham Yearwood, PhD¹, Gina DeStefano, MS¹, Helen Koenig, MD,MPH^{2,3}

1. UrSure Inc 2. Philadelphia FIGHT 3. Perelman School of Medicine at the University of Pennsylvania



Background

- There are 1.2 million people at high-risk of HIV infection in the US and millions worldwide.
- Tenofovir (TFV) is the major component of Pre-Exposure Prophylaxis (PrEP) for HIV prevention and a critical component of antiretroviral therapy (ART) for the treatment of HIV.
- For both the prevention and treatment of HIV, adherence is critical.
- Current adherence to TFV-containing regimens is sub-optimal, but effective methods for drug monitoring can improve adherence.
- There is a need for an affordable, disposable, easy to use point of care (POC) test to assess adherence in order to better manage PrEP and ART care.

Methods

Synthesized a novel TFV derivative to be conjugated to KLH and be used as an immunogen to produce anti – TFV monoclonal antibodies

Screened and selected antibodies in a competitive enzyme linked immunoassay (ELISA) format to assess sensitivity and specificity

LFIA prototype was developed on a nitrocellulose pad with a sample port, antibody gold conjugate strip, immobilized TFV-BSA strip and a control


Varying spiked TFV concentration samples were added to the prototype and sensitivity was assessed

Objective

To establish a qualitative lateral flow immunoassay (LFIA) to measure the concentration of TFV in a patient's urine to measure and promote adherence.

Lateral Flow Immunoassay

Strip Components




- The sample migrates upward starting from the sample port and ending at the control line.
- First, the sample is added to the port and then migrates to the antibody gold conjugate.
- If TFV is present in the sample, it will bind to the conjugate and not be captured at the immobilized TFV-BSA conjugate test line = no 1st line.
- If TFV is not present in the sample, the antibody conjugate will not bind and it will be captured at the immobilized TFV-BSA conjugate test line = 1st line.
- The control line (C) is the built-in procedural control, which consists of an anti-mouse IgG = 2nd line.

Results

mAB Sensitivity and Specificity of ELISA

	LC-MS (+)	LC-MS (-)
Antibody (+)	50	8
Antibody (-)	0	140

LFIA Prototype Standard Curve



- The monoclonal antibody performance in the competitive ELISA format demonstrated 100% sensitivity and 94.6% specificity.
- When the TFV-BSA strip was dipped in the gold nanoparticle antibody conjugate the results showed sensitive detection of TFV spiked samples across a dilution series.
- The prototype has a limit of detection of 1,250 ng/mL.

Conclusion

A POC LFIA prototype that could detect TFV concentrations in the urine matrix was developed and could prove to be a useful tool in promoting optimal adherence for HIV prevention and treatment.

Future Steps

- Optimization, validation, verification, and manufacturing initiation of the POC LFIA test.
- Clinical trials to assess the acceptability, feasibility, and effectiveness of the test in real world settings.
- Scale-up of the POC test to patients and providers worldwide.

References

1. HIV Surveillance Report, 2014 [Internet]. Centers for Disease Control and Prevention.; 2016 Nov [cited 2017 Aug 12]. Available from: <http://www.cdc.gov/hiv/library/reports/surveillance/>. Published November 2015
2. FACT SHEET JULY 2017 [Internet]. UNAIDS; 2017 Jul. Available from: http://www.unaids.org/sites/default/files/media_asset/UNAIDS_FactSheet_en.pdf
3. Number of Deaths Due to HIV/AIDS. [Internet]. World Health Organization Global Health Observatory; 2016 Apr. Available from: http://www.who.int/gho/hiv/epidemic_status/deaths/en/
4. Kearney BP, Flaherty JF, Shah J. Tenofovir disoproxil fumarate: clinical pharmacology and pharmacokinetics. Clin Pharmacokinet. 2004;43(9):595–612.
5. Amico KR, Marcus JL, McMahan V, Liu A, Koester KA, Goicochea P, et al. Study product adherence measurement in the iPrEx placebo-controlled trial: concordance with drug detection. J Acquir Immune Defic Syndr 1999. 2014 Aug 15;66(5):530–7.
6. Golin CE, Liu H, Hays RD, Miller LG, Beck CK, Ickovics J, et al. A prospective study of predictors of adherence to combination antiretroviral medication. J Gen Intern Med. 2002 Oct;17(10):756–65.